

EPROUVETTES CLASSE B

Réf	Description	volume ml	divisionl ml	tolerance ml	couleur graduation	Qté / carton	Prix HT € / carton
2230154	5 : 0,1 ml	5	0,1	0.075	bleue	2	11,00
2230160	10 : 0,2 ml	10	0,2	0.15	bleue	2	11,00
2230170	25 : 0,5 ml	25	0,5	0.375	bleue	2	12,00
2230175	50 : 1 ml	50	1	0.75	bleue	2	12,25
2230180	100 : 1 ml	100	1	0.75	bleue	2	13,70
2230186	250 : 2 ml	250	2	1.5	bleue	2	30,00
2230190	500 : 5 ml	500	5	3.75	bleue	2	42,00
2230193	1000 : 10 ml	1000	10	7.5	bleue	1	33,20
2230194	2000 : 20 ml	2000	20	15	bleue	1	57,40




EPROUVETTES CLASSE A

Réf	Description	volume ml	divisionl ml	tolerance ml	couleur graduation	Qté / carton	Prix HT € / carton
2240154	5 : 0,1 ml	5	0,1	0.05	bleue	2	14,10
2240160	10 : 0,2 ml	10	0,2	0.1	bleue	2	13,90
2240170	25 : 0,5 ml	25	0,5	0.25	bleue	2	14,80
2240175	50 : 1 ml	50	1	0.5	bleue	2	16,00
2240180	100 : 1 ml	100	1	0.5	bleue	2	17,40
2240186	250 : 2 ml	250	2	1	bleue	2	32,70
2240190	500 : 5 ml	500	5	2.5	bleue	2	56,30
2240193	1000 : 10 ml	1000	10	5	bleue	1	41,20
2240194	2000 : 20 ml	2000	20	10	bleue	1	79,60



**FIOLES JAUGEES CLASSE A VERRE DURAN
AVEC BOUCHON POLYETHYLENE**

Réf	Description	volume ml	tolérance ml	hauteur mm	couleur graduation	Qté / carton	Prix HT € / carton
2820152	5 ml NS 7/16	5	0.025	70	bleue	2	23,00
2820153	5 ml NS 10/19	5	0.04	70	bleue	2	23,50
2820160	10 ml NS 7/19	10	0.025	90	bleue	2	23,00
2820161	10 ml NS 10/19	10	0.04	90	bleue	2	23,60
2820165	20 ml NS 10/19	20	0.04	110	bleue	2	24,50
2820166	20 ml NS 12/21	20	0.06	105	bleue	2	24,50
2820170	25 ml NS 10/19	25	0.04	110	bleue	2	24,50
2820171	25 ml NS 12/21	25	0.06	110	bleue	2	27,00
2820175	50 ml NS 12/21	50	0.06	140	bleue	2	28,00
2820178	50 ml NS 14/23	50	0.08	140	bleue	2	28,00
2820180	100 ml NS 12/21	100	0.1	170	bleue	2	29,50
2820181	100 ml NS 14/23	100	0.1	170	bleue	2	29,50
2820185	200 ml NS 14/23	200	0.15	210	bleue	2	34,00
2820186	250 ml NS 14/23	250	0.15	220	bleue	2	36,00
2820190	500 ml NS 19/26	500	0.25	260	bleue	2	41,60
2820192	1000 ml NS 29/3	1000	6	300	bleue	1	36,30
2820193	1000 ml NS 24/2	1000	0.4	300	bleue	1	33,45
2820194	2000 ml NS 29/3	2000	0.6	370	bleue	1	45,50
2820197	5000 ml NS 34/3	5000	1.2	470	bleue	1	94,60
2820198	10000 ml NS 45/	10000	2	570	bleue	1	182,00



**FIOLES JAUGÉES CLASSE A VERRE DURAN
AVEC BORD VERSEUR**

Réf	Description	volume ml	tolérance ml	hauteur mm	couleur graduation	Qté / carton	Prix HT € / carton
2800152	5 ml	5	0.025	70	blue	2	15,00
2800160	10 ml	10	0.025	90	blue	2	15,00
2800165	20 ml	20	0.04	110	blue	2	15,00
2800170	25 ml	25	0.04	110	blue	2	15,00
2800175	50 ml	50	0.06	140	blue	2	17,40
2800180	100 ml	100	0.1	170	blue	2	20,00
2800185	200 ml	200	0.15	210	blue	2	27,00
2800186	250 ml	250	0.15	220	blue	2	28,40
2800190	500 ml	500	0.25	260	blue	2	36,60
2800193	1000 ml	1000	0.4	300	blau	1	26,05
2800194	2000 ml	2000	0.6	370	blau	1	40,75
2800197	5000 ml	5000	1.2	470	blue	1	88,80
2800198	10000 ml	10000	2	570	blue	1	172,00